

## PH30P Solid Piercing pH Meter For Meat, fruits, soil, etc..



- Grab handle shell for your stable and comfortable grip.
- pH Range: 0-13pH
- Strong and sensitive pH bulb, accurate and durable.
- Removable 316L stainless steel blade, easy cleaning.
- Waterproof and dustproof housing, IP67 rated.
- Precision easy operation: all functions operated in one hand.
- Large LCD with backlight.
- Real time electrode efficiency icon indication.
- 1\*1.5 AAA long battery life.
- Auto-Power Off saves battery after 5 mins non-use.
- Auto Lock Function

Range	0.00 ~ +13.00 pH
Resolution	0.01 pH
Accuracy	± 0.01 + 1LSD
Temp. Range	0 - 100.0 °C / 32 - 212 °F
Working Temp.	0 - 60.0 °C / 32 - 140 °F
Temp. Resolution	0.1 °C / 1 °F
Calibration Point	Up to 3 (auto recognized)
pH Buffer Kinds	USA - pH 4.01, 7.00, 10.01 NIST - pH 4.01, 6.86, 9.18
Electrode	CS1930 Solid Piercing pH electrode (316L Stainless Piercing Blade)
LCD	20*30 mm multi-line crystal display with backlight
Piercing depth	≥ 28mm
Temp. Compensation	Auto / Manual
Lock	Auto / Manual
Auto backlight off	1 minute
Auto power off	5 minutes without key be pressed
Power	1x1.5V AAA battery
Dimensions	(H×W×D) 185×40×48 mm
Weight	95g
Protection	IP67

